

DataTrace Pro Database

If your DataTrace Pro Database is getting too large, you may be receiving the following error message when you are trying to view a study, even when study names are present:

“No Data Found Expand Time Range.”

This means that your database may be full. You should archive data using the archive utility within DT Pro until the database reaches a manageable size.

If you encounter any errors, try archiving a smaller set of data at a time.

You will also want to check the following:

- DT Pro Application Data Folder
 - DTData.mdf should be less than 4GB
 - To access the Folder in Vista through Windows 8, go to:
C:\ProgramData\DTPro\Data
 - To access the Folder in CP, go to:
C:\Documents and Settings\All Users\Application Data\DTPro\Data

