

DTRF Release Notes:

Current Version: 1.1.8; Build Date: 10/12/2009

This document contains historical and important information about the release of DTRF packages.

RELEASE HISTORY

=====

Rev 1.1.5 Initial Release

Rev 1.1.6 Bug Fixes

Corrected bug pertaining to saving user settings when security was enabled
Corrected oversight of User message in Test Tracer to identify detection of

Low Battery Radio Reset

Adjusted maximum allowed radio interval to 15 minutes

Disabled menu item for unsupported thermal map printing functionality

Correct error handling to allow diagnostic reporting for tracer's with no

data

Correct oversight which made invisible the Read Tracers Batch View > History

Tab > Row Omission menu item

Rev 1.1.7 Bug Fixes

Corrected communications problems when software is used on an Asian computer
(2 byte text memory format)

separator for other countries

Corrected several issues pertaining to date/time display format and decimal

signature per Read session (match DTW 4.05 functionality)

Corrected bug pertaining to RF data recovery after logger memory is full

Enabled RF Data Encryption Mode

Corrected battery display > 100% issue

Corrected max. readings for MPRF type in diagnostic report

old-style DB9 COM ports

Corrected bug pertaining to newer laptops or computers that do not have the

digits

Display format for pressure display settings can now be set for 3 decimal

tracer intervals

Set New Battery functionality is now available for MPIII loggers

Corrected Stop On Time issue so that Stop On Time will now work for all

corruption

Corrected the corruption of MPRF Stop Bytes that resulted in a header

Rev 1.1.8 Bug Fixes

Enabled the ability to turn off all MPRF related functions

security level has different areas of access

Modified code to support security functionalities similar to DTW, i.e. each

locked out

Resolved issue where logger reading displays in Placement view were not

database

Resolved issue where custom lethality settings were not saved to the

data for the Host PCIF

Resolved issue where the Programming Report screen was reporting incorrect

Fixed issues with User Calibration for both MPRF and MPIII loggers

Added a message to warn users when they are about to close Batch View screen

Resolved issue with User Comment being persistent when a logger is

programmed

Resolved issue with Run ID Length option box should only be visible for

Double-byte OSes

Resolved issue with new MPIII single channel temp. loggers that can store

more than 16000 data points

Fixed issue with writing tracer start time to the audit trail DB when a

DTRF Release Notes.txt

tracer is programmed

Fixed issue with Batch Read mode timed out when reading FRB tracers

Fixed issue with not showing correct Run IDs on the Program List, which is the list of all tracers that have just been programmed

Fixed issue to add an entry to the Audit Trail DB when an out-of-range error is detected upon reading a tracer

Added new code to make sure that when decimal settings changed in DTRF are reflected in Mini-DTW so reports can pick them up

Resolved issue with current logged in username does not show up on a printed report

Fixed issue with Calibration dates intermittently showing up in tracer diagnostic reports

Added new feature to support the cropping of %RH values on display (from 0 to 100)

INSTALLER / CD MODIFICATIONS

=====

Rev 1.1.6

Added Copy of Manual to install directory

Added this release notes file to CD, with copy to install directory

Rev 1.1.7

Suppressed installation messages relating to existence of O.S. shared DLLs.

Included the latest digitally-signed version of USB driver.

Added upgrade instructions PDF to CD.

VALIDATION INFORMATION

=====

Rev 1.1.5 DTRF Validation Master v1.0, Full implementation, 6/16/08

Rev 1.1.6 DTRF Validation Master v1.1, Partial implementation, 08/11/08

Section 8.1, with Security enabled, verify Temp units storage.

Section 4.3.13 and 4.3.14, check 15 minute maximum interval.

Section 3.4.1, using a logger with no data

Verified Row omission is available

Verified Thermal Map print menu is not available.

Rev 1.1.7 DTRF Validation Master v1.2, Full implementation, 11/1/08

KNOWN ISSUES

=====

1. Batch View functionality in Read window is not documented in Help or the User Manual

2. User message to identify detection of Low Battery Radio Reset is not documented in Help or the User Manual and is English only.

3. Software will not function properly on a Vista O.S. PC when installed by one logged in user and used by a different logged in user.

4. In Full security mode reporting does not include automatic annotation of logged in user on reports

5. Alarm settings currently does not support pressure tracers. In addition, RH alarm limits are currently being used for all types of second channel data.