



MesaLabs

DataTrace Customer,

DataTrace MPIII Dataloggers have been tested and used by the industry for several years for many applications including validations. DataTrace is suitable for temperature mapping and validation in steam sterilizers, ovens and other appliances over the calibrated range of -40C to 140C. DataTrace Temperature MPIII's in particular have a response time that is appropriate for the standard logging interval of 1 to 2 minutes. (Please see the attached document for further response time information.

DataTrace Product Manager

DataTrace Response to Change

The speed at which a sensor responds to a change in temperature can be critical to monitoring a process accurately. Every temperature sensor requires a certain time for the temperature reported to react to, and equilibrate with, a change in ambient temperature. Two terms are frequently used to describe this analysis: Thermal Time Constant and Response Time. They are frequently used interchangeably to describe this analysis; however, they are not the same.

The following is the technical description according to MIL-T-23648A:

Thermal Time Constant: The Thermal Time Constant is the time required for a temperature sensor to change 63.2% of the total difference between its initial temperature and final body temperature when subjected to a step function change in temperature.

As defined above, the thermal time constant, t_{th} , is the time that a sensor takes to react to a rapid change in the ambient temperature and to reach 63.2% of the actual temperature change. Typically, in order to determine this value, a sensor is moved rapidly from one temperature environment to another. The time taken to reach the second value, or Thermal Time Constant, is then measured.

In many applications, the Thermal Time Constant is not the most important measure of response. More important for many users of temperature sensors is Response Time.




Response Time, or t_{90} , is the time required for a temperature sensor to achieve 90% of the actual temperature change.

Typically, t_{90} is slightly more than two times the t_{th} value for the same sensor configuration in the same environment.

The actual reaction time depends on the size of the sensor, the media in which the temperature is being measured, and on how good the thermal contact is between the sensor and the target. Small sensors in fluid react quickest and large sensors in air react slowest. Sensors react 10 to 50 times faster in a well-stirred fluid than in still air.

Both the Thermal Time Constant and the Response Time are usually more important in a rapidly changing process environment than in a more stable storage environment. However, they are important concepts to be aware of as measurement needs are considered.

DataTrace Technical Note

PROBE	SENSOR	PROBE TIP DIAMETER	WELL-STIRRED FLUID		AIR
			63%	90%	63%
Rigid Bullet Tip 	Thermistor	0.188" (4.78mm)	8 sec.	21 sec.	36 sec.
Stainless Steel Taper Tip (EXT) 	RTD	0.125" (3.18mm)	<8 sec.	17 sec.	26 sec.
Bendable Flex Rounded Tip 	RTD	0.079" (2.00mm)	<3 sec.	<6 sec.	24 sec.
Rigid Taper Tip Stainless Steel Flex Taper Tip Silicone Flex Taper Tip 	Thermistor	0.125" (3.18mm)	2 sec.	5 sec.	27 sec.
Stainless Steel Flex Needle Tip Silicone Flex Needle Tip 	Thermistor	0.058" (1.5mm)	<2 sec.	<5 sec.	11 sec.