



MesaLabs

DataTrace[®] Calibration Warranties (US Only)

Mesa Laboratories, Inc. provides a warranty on our NIST traceable calibration services, in addition to our 12 month hardware warranty on new products. To avoid confusion on the Calibration Warranties for DataTrace Loggers we are including the following for your information. Note that factory calibration warranty can vary depending on type of Logger and whether the subject Logger is a new build or a Logger being returned from service.

- All NEW Temperature and Pressure Loggers and the Temperature component of Humidity Loggers are warranted to hold specification for 12 months from date of original shipment.
- All SERVICED Temperature and Pressure Loggers and the Temperature component of Humidity Loggers are warranted to hold specification for 6 months.
- All NEW and SERVICED Humidity Loggers are warranted to hold the Humidity component specification for 3 months.

Field calibration of the humidity sensor is possible on all Humidity Loggers. The Calibration Utility is built into the DataTrace Pro software program and the instructions for this function are available in the help section under “How to Calibrate loggers or monitors”.