

# PC Interface

# DataTrace MPIII/MPRF PC Interface

The DataTrace PC Interface is an integral part of the DataTrace data logger system. Mesa's USB compatible PC interface rapidly communicates with your data loggers using wireless infrared (IR) technology. Communication takes only seconds, and using wireless 2.4GHz technology, the MPRF interface is also capable of communicating with large groups of MPRF data loggers in real-time, connecting DataTrace data loggers with Mesa's powerful DT Pro software.



MPIII

**Specifications:**

<b>Height:</b>	1.44"
<b>Length:</b>	4.30"
<b>Width:</b>	2.75"

**Specifications:**

<b>Height:</b>	1.44"
<b>Length:</b>	4.30"
<b>Antenna Length: (MPRF)</b>	4.42"
<b>Width:</b>	2.75"



MPRF