

## Helpful MPIII Tips

DataTrace Tech Support:  
DTSupport@MesaLabs.com



New to DataTrace? Here are some things you might like to know!

### Get to know your logger! Temperature parameters are included.

- Temperature: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPIII-Temperature-Brochure.12DEC2017.pdf>
- Pressure: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPIII-Pressure-Brochure.20JAN2016.pdf>
- Humidity: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPIII-Humidity-Brochure.23NOV2015.pdf>

### When it comes time for logger calibration:

Please fill out the form linked below and email it to us. We'll provide a quote from there.

- <http://datatrace.mesalabs.com/calibration/>
- [CustomerService@MesaLabs.com](mailto:CustomerService@MesaLabs.com)

### Need to purchase new equipment? Your regional representative will be happy to help with quotes.

- [DataTraceSales@MesaLabs.com](mailto:DataTraceSales@MesaLabs.com)

### Working at high temperatures? Logger Armor will help protect your logger and prolong its life:

- <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/DT-Thermal-Packs-Brochure.08APR2015.pdf>

### We offer an MPIII logger wrench for easy opening and closing of the battery cap. This will help protect against fluid intrusion:



MPIII Battery Cap Wrench  
Part number DT170002

## **Battery Changes:**

- A PDF of battery change instructions can be found here: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2015/03/MPIII-Battery-Change-Quick-Start-Guide.05OCT2016.pdf>
- A video is here: <https://www.youtube.com/watch?v=FBk9nIFO8N8>
- Please note that it is important to change the o-ring every time the battery cap is opened. The o-ring acts as a sacrificial seal and protects against damaging fluid intrusion.
- New MPIII batteries come in packs of 20 with silicone grease and o-rings. Part number 124804-500.

## **Software license information:**

- Each software license is good to activate the software on one computer.
- Please let us know if you need to switch the software to a new computer and we will be happy to provide detailed instructions.

**Do you need help with your DT Pro Software? Navigate to the Help menu at the top of your screen where you'll find a searchable manual as well as a step-by-step tutorial.**

**Pressing F1 in ANY screen of DT Pro will bring up the applicable Help topic.**