

Helpful MPRF Tips

DataTrace Tech Support:
DTSupport@MesaLabs.com



New to DataTrace? Here are some things you might like to know!

Get to know your logger! Temperature parameters are included.

- Temperature: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPRF-Temperature-Brochure.12MAY2017.pdf>
- Pressure: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPRF-Pressure-Brochure.12MAY2017.pdf>
- Humidity: <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/MPRF-Humidity-Brochure.24JUL2017.pdf>

When it comes time for logger calibration:

- Please fill out the form linked below and email it to us. We'll provide a quote from there.
<http://datatrace.mesalabs.com/calibration/>
CustomerService@MesaLabs.com

Need to purchase new equipment? Your regional representative will be happy to help with quotes.

- DataTraceSales@MesaLabs.com

We offer an MPRF logger wrench for easy opening and closing of the battery cap. This will help protect against fluid intrusion:



MPRF Battery Cap Wrench
Part number 270009-1

Working at high temperatures? Logger Armor will help protect your logger and prolong its life:

- <http://datatrace.mesalabs.com/wp-content/uploads/sites/19/2014/01/DT-Thermal-Packs-Brochure.08APR2015.pdf>



Logger Armor

Working at cryogenic temperatures? Special batteries will be needed if you're working with low-temperature loggers.

- <http://datatrace.mesalabs.com/2014/06/06/80-low-temp-batteries-for-dataloggers/>

Changing a Battery:

- Battery change instructions can be found here: <https://datatrace.mesalabs.com/mprf-battery-change-quick-start-guide/>
- A video is here: https://www.youtube.com/watch?v=NvQ41-R_ZaQ
 - Please note that it is important to change the o-ring every time the battery cap is opened. The o-ring acts as a sacrificial seal and protects against damaging fluid intrusion.
- Part numbers and other purchasing details are here: <http://datatrace.mesalabs.com/datatrace-mprf-accessories/>

Software license information:

- Each software license may be activated on one computer.
- Please let us know if you need to switch the software to a new computer and we will be happy to provide detailed instructions.

Do you need help with your DT Pro Software? Navigate to the Help menu at the top of your screen where you'll find a searchable manual as well as a step-by-step tutorial.

Pressing F1 in ANY screen of DT Pro will bring up the applicable Help topic.